

FuelCell Addendum – Arbitrary Control Termination

D. Johnson, Scribner Associates, Inc.

4/16/2007, Ver. 2

Introduction

This addendum describes the use of new termination actions in the Arbitrary Control Experiment in FuelCell Version 3.8d and later.

The termination is triggered only when data points are recorded using Action 0.

The experiment termination is controlled through Actions 101-116

```
' Action = 101      terminate when Current < value
' example         101    0.5    7
' terminate if the current is less than 0.5 amps
' first parameter is the Current Limit in Amps
' second parameter controls the type of termination
' 0 = no termination
' 1 = terminate the Arbitrary Control Experiment, continues with the next experiment in
the list
' 2 = does not perform experiments after the Arbitrary Control Experiment and before the
end of a repeat loop
' 4 = exits a repeat loop, continues with experiments after the end of the loop
' 8 = stops all experiments in the experiment list
' The termination types are additive: 7 = 1+2+4 = terminate Arbitrary Control Experiment
and all experiments until the end of the loop. Exit the loop and perform experiments
after the end if the loop
'
' Action = 102      terminate when Current > value
' example         102    0.5    7
' terminate if the current is greater than 0.5 amps
' first parameter is the Current Limit in Amps
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 103      terminate when Potential < value
' example         103    0.5    7
' terminate if the potential is less than 0.5 volts
' first parameter is the Voltage Limit in Volts
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 104      terminate when Potential > value
' example         104    0.5    7
' terminate if the potential is greater than 0.5 volts
' first parameter is the Potential Limit in Volts
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 105      terminate when Power < value
' example         105    0.5    7
' terminate if the power is less than 0.5 watts
' first parameter is the Power Limit in Watts
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 106      terminate when Power > value
' example         106    0.5    7
' terminate if the power is greater than 0.5 watts
```

```

' first parameter is the Power Limit in Watts
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 107          terminate when iR Corrected Potential < value
' example              107    0.5    7
' terminate if the potential is less than 0.5 volts
' first parameter is the Voltage Limit in Volts
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 108          terminate when iR Corrected Potential > value
' example              108    0.5    7
' terminate if the potential is greater than 0.5 volts
' first parameter is the Potential Limit in Volts
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 111          terminate when Cell Temperature < value
' example              105    65     7
' terminate if the temperature is less than 65 C
' first parameter is the Temperature Limit in degrees C
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 112          terminate when Cell Temperature > value
' example              106    0.5    7
' terminate if the temperature is greater than 65 C
' first parameter is the Temperature Limit in degrees C
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 113          terminate when Anode Temperature < value
' example              105    65     7
' terminate if the temperature is less than 65 C
' first parameter is the Temperature Limit in degrees C
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 114          terminate when Anode Temperature > value
' example              106    0.5    7
' terminate if the temperature is greater than 65 C
' first parameter is the Temperature Limit in degrees C
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 115          terminate when Cathode Temperature < value
' example              105    65     7
' terminate if the temperature is less than 65 C
' first parameter is the Temperature Limit in degrees C
' second parameter controls the type of termination
' see Action 101 for a description of the termination types
'
' Action = 116          terminate when Cathode Temperature > value
' example              106    0.5    7
' terminate if the temperature is greater than 65 C
' first parameter is the Temperature Limit in degrees C
' second parameter controls the type of termination
' see Action 101 for a description of the termination types

```

Termination Types and Examples

Several different termination types are available. Each type terminates or skips different parts of the experiment list. They can be combined to create more complex terminations.

In the examples, the voltage at 2 Amps is less than 0.5 V and triggers the termination

Termination Type 1 Example:

Skips the remaining lines in the Arbitrary Control Experiment. Other experiments in the list, and the Repeat Loop are not skipped.

```

Begin Loop
  Arbitrary Control Experiment
    103  0.5  1  \terminate if potential < 0.5, termination type 1
    0     5      \take point
    1     1      \set current to 1 Amp
    0     5      \take point
    1     2      \set current to 2 Amp
    0     5      \take point
    1     1      \set current to 3 Amp
    0     5      \take point
  Open Circuit Experiment
  Potentiostatic Experiment
End Loop
Begin Loop
  Arbitrary Control Experiment
  Open Circuit Experiment
  Potentiostatic Experiment
End Loop

```

Termination Trigger (points to the line: \set current to 2 Amp)

Skip remaining lines in the Arbitrary Control Experiment (points to the line: \set current to 3 Amp)

Termination Type 2 Example:

Skips the remaining experiments until the end of the Repeat Loop. The remaining lines in the Arbitrary Control Experiment are performed. The Repeat Loop will continue to operate.

```

Begin Loop
  Arbitrary Control Experiment
    103  0.5  2  \terminate if potential < 0.5, termination type 2
    0     5      \take point
    1     1      \set current to 1 Amp
    0     5      \take point
    1     2      \set current to 2 Amp
    0     5      \take point
    1     1      \set current to 3 Amp
    0     5      \take point
  Open Circuit Experiment
  Potentiostatic Experiment
End Loop
Begin Loop
  Arbitrary Control Experiment
  Open Circuit Experiment
  Potentiostatic Experiment
End Loop

```

Termination Trigger (points to the line: \set current to 2 Amp)

Skip remaining experiments until the end of the loop (points to the line: Potentiostatic Experiment)

Termination Type 4 Example:

Exits the Repeat Loop. The remaining lines in the Arbitrary Control Experiment and other experiments in the list, and the Repeat Loop are not skipped.

```

Begin Loop
  Arbitrary Control Experiment
    103  0.5  4  \terminate if potential < 0.5, termination type 4
    0    5    \take point
    1    1    \set current to 1 Amp
    0    5    \take point
    1    2    \set current to 2 Amp
    0    5    \take point
    1    1    \set current to 3 Amp
    0    5    \take point
  Open Circuit Experiment
  Potentiostatic Experiment
End Loop
Begin Loop
  Arbitrary Control Experiment
  Open Circuit Experiment
  Potentiostatic Experiment
End Loop

```

Termination Trigger

Exit the Repeat Loop

Termination Type 8 Example:

Ends all Experiments. The remaining lines in the Arbitrary Control Experiment are not skipped.

```

Begin Loop
  Arbitrary Control Experiment
    103  0.5  8  \terminate if potential < 0.5, termination type 8
    0    5    \take point
    1    1    \set current to 1 Amp
    0    5    \take point
    1    2    \set current to 2 Amp
    0    5    \take point
    1    1    \set current to 3 Amp
    0    5    \take point
  Open Circuit Experiment
  Potentiostatic Experiment
End Loop
Begin Loop
  Arbitrary Control Experiment
  Open Circuit Experiment
  Potentiostatic Experiment
End Loop

```

Termination Trigger

Skip all remaining experiments in the Experiment List

Combining Termination Types:

Several Termination types may be combined.

Termination Type 7 Example:

Termination Type 7 combines the results of types 1, 2 and 4 (1+2+4=7)

This Type will skip all remaining lines in the Arbitrary Control Experiment, skip all experiments between the Arbitrary Control Experiment and End Loop, and exit the Repeat Loop.

```

Begin Loop
  Arbitrary Control Experiment
    103  0.5  7  \terminate if potential < 0.5, termination type 7
    0    5          \take point
    1    1          \set current to 1 Amp
    0    5          \take point
    1    2          \set current to 2 Amp
    0    5          \take point
    1    1          \set current to 3 Amp
    0    5          \take point
  Open Circuit Experiment
  Potentiostatic Experiment
End Loop
Begin Loop
  Arbitrary Control Experiment
  Open Circuit Experiment
  Potentiostatic Experiment
End Loop

```

Termination Trigger

Skip all remaining lines in the Arbitrary Control Experiment, skip remaining experiments in the Loop, exit Repeat Loop

Termination Type 0:

Termination Type 0 can be used to turn off the termination condition.

For example, in the following Arbitrary Control list, the termination condition will only be checked during the data point measured at 2 Amps. The other data points will not trigger a termination.

```

Arbitrary Control Experiment
  0    5          \take point
  1    1          \set current to 1 Amp
  0    5          \take point
  1    2          \set current to 2 Amp
  103  0.5  7    \terminate if voltage < 0.5
  0    5          \take point
  103  0.5  0    \turn off termination
  1    3          \set current to 3 Amp
  0    5          \take point

```