



Membrane Characterization in Electrolyte

K.R. Cooper, Scribner Associates Inc., Z. Wang, J. Parrondo, V. Ramani,

MOTIVATION

- New ionomer membranes and separators are being developed for emerging electrochemical processes especially in the field of energy storage and conversion such as **REDOX FLOW BATTERIES**
- A commercial cell for *ex situ* membrane characterization in relevant environment is needed

OBJECTIVE

Develop a versatile tool to support membrane R&D used with liquid environments. Desired features:

- 4-electrode to eliminate non-ohmic impedance sources
- Small electrolyte volume
- Small sample size
- Easy assembly/disassembly
- Compatible with wide range of environments

SOLUTION: Commercial H-Cell

Design Features

- 4-electrode with removable REs
- High-surface area Pt CE & WE
- 30 mL electrolyte volume per half
- Ports for gas inlet/outlet & Temp
- Adjustable Luggin probe position
- Chemically resistant materials

Correcting $R_{measured}$ for electrolyte resistance, $R_{electrolyte}$


$$\sigma_{membrane} (S/cm) = \frac{t}{R_{membrane} \cdot A}$$

t = membrane thickness (cm)
 A = membrane area (cm²)

- Resistance R (Ω) determined using LSV (± 0.1 V vs. OCV, 50 mV/s)
- $R_{measured} = R_{membrane} + R_{electrolyte}$
- $R_{electrolyte}$ = Uncompensated electrolyte resistance
- $R_{electrolyte}$ obtained by
 - Measuring H-cell resistance without a membrane
 - Calculating from electrolyte conductivity $\sigma_{electrolyte}$ (S/cm) and the cell constant K_{cell} (cm⁻¹)


$$R_{electrolyte} = \frac{K_{cell}}{\sigma_{electrolyte}}$$

- K_{cell} depends only on cell geometry (Luggin probe location, cell cross-section) and is independent of $\sigma_{electrolyte}$
- Conductivity standards used to obtain K_{cell}




EVALUATION OF NOVEL AEM FOR RFBs

- Polymer backbone: polystyrene-*block*-poly-(ethylene-*ran*-butylene)-*block*-polystyrene (SEBS)
- Functional group: trimethylamine (TMA)




SEBS-TMA Properties

| | |
|---------------------------------------|--------------|
| IEC (mmol/g) | 1.35 ± 0.02 |
| Water uptake | 52% |
| Permselectivity | 73% |
| Transport no. ($t_{Cl^-}; t_{K^+}$) | (0.87; 0.13) |
| Tensile strength (MPa) | 1.83 |
| Elongation at break | 350% |




- AEM σ in Cl^- > in SO_4^{2-}
- σ increases slightly over time



- AEM σ in Cl^- > in SO_4^{2-}
- $\sigma_{SEBS-TMA} \sim 1/3 \sigma_{Nafion}$

- σ : TP (H-cell in Cl^-) > IP (BT-110 in high purity H₂O)



CONCLUSIONS

- *H-cell* is a viable cell design for reproducible through-thickness resistance and conductivity measurement of ion exchange membranes and separators in liquid electrolytes
- 4-electrode measurement allows simple DC measurement (LSV) with straight forward data analysis of resistance
- Reproducible, accurate results obtained with Nafion and an experimental AEM in multiple RFB-type electrolytes