

Characterizing Through-Plane & In-Plane Ionic Conductivity of Polymer Electrolyte Membranes

Kevin R. Cooper
Scribner Associates, Inc.

Paper # B10-1043 presented at PEMFC-11, *220th ECS Meeting, Boston, MA*

October 12, 2011

Portions of this work were supported by the U.S. DOE through the University of Central Florida (# DE-FC36-06GO16028). Data provided by BekkTech LLC is greatly appreciated.


Orientation and Nomenclature

TP = through-plane or through-thickness

IP = in-plane

L = longitudinal = parallel to machine (extrusion) direction



T = transverse = perpendicular to machine direction



Is the conductivity of Nafion® isotropic? ... No consensus in published literature

- No, it is anisotropic
 - $\sigma_{IP}:\sigma_{TP} = 3.6$ [Gardner]
 - $\sigma_{IP}:\sigma_{TP} = 2.5 - 5$ (w/ Pressure) [Ma]
 - $\sigma_{IP}:\sigma_{TP} = 1.8 - 5$ [Casciola]
 - ✓ This work (N112)
- Yes, it is isotropic, $\sigma_{IP}:\sigma_{TP} \approx 1$
 - Nouel, Silva
 - ✓ Cooper (NR-212)
- Discrepancy due
 - Different water content (λ)
 - Not accounting for non-membrane ohmic contributions in TP measurements
- Anisotropic within the plane for extruded material: $\sigma_{IP-Long} > \sigma_{IP-Trans}$


Casciola *et. al.*, *J. Power Sources* **162** 141 (2006)
 Cooper, *JES* **157**, B1731 (2010)
 Gardner *et. al.*, *J. Electroanal Chem* **449** 209 (1998)
 Ma *et. al.*, *JES* **153** A2274 (2006)
 Nouel, Fedkiw, *Electrochimica Acta* **43** 2381 (1998)
 Silva, *et. al.*, *J. Power Sources* **134** 18 (2004)



Objective: Measure and compare in-plane and through-plane ion conductivity of polymer electrolyte membranes

➤ In-plane

- Pro: geometry allows true 4-electrode configuration → implementation & data interpretation are easy
- Con: not orientation of interest




$$\sigma_{IP} = \frac{L}{R \cdot A} = \frac{L}{R \cdot W \cdot T}$$

➤ Through-plane

- Pro: measured parameter is in orientation of interest
- Con: measurement and data analysis more difficult
- Instrument & methods recently developed [1]

$$\sigma_{TP} = \frac{T}{R \cdot A}$$



1. Cooper, *J. Electrochem. Soc.* **157**, B1731 (2010)

Experimental

- Materials
 - Nafion[®] 112 and NR-212 – non-supported PEM
 - GORE-SELECT[®] – supported PEM
- Through-plane – Membrane Test System *MTS 740* [1]
 - Integrated membrane clamp and electrodes
 - Environmental control & measurement – T, dew point, RH, gas flow, *etc.*
- In-plane – BekkTech BT-110 Conductivity Clamp
 - Located in MTS 740 chamber for environmental control

1. Cooper, *J. Electrochem. Soc.* **157**, B1731 (2010)



GORE-SELECT, GORE and designs are trademarks of W. L. Gore & Associates, Inc.

Procedure – Pre-test

- As-received membrane, stored at ambient conditions
- 32 mm x 10 mm sample
- Measure “dry” membrane thickness
 - Mean of 5 locations, 3x measurements/location
 - Brunswick Instruments
 - low load (50 g)
 - high accuracy gage ($\pm 0.2 \mu\text{m}$ for $50 \mu\text{m}$ sample)








Brunswick Instruments Film Thickness Gage

Procedure


- Temperature series (°C): 80 → 30 → 120
- Per temperature
 - Wet-up 2 hr @ 70% RH
 - RH cycle: 70 → 20 → 95 %, 15 min or 95% → 20%, 30 min
 - Resistance measurement after 15 or 30 min
- Through-plane: 2-electrode / 4-terminal impedance measurement
 - 10 MHz – 1 Hz, 10 mV_{AC}, 0 V_{DC}
- In-plane: 4-electrode / 4-terminal DC measurement

Temp, °C	Total Dry Gas Flow, sccm	Pressure, kPa _a
30	500	100
80	500	100
120	500	230

Post-test Procedure – EIS Analysis of Through-Plane Data



Through-thickness Measurement Method and Ohmic Resistance Sources



➤ Ohmic resistances that contribute to the high frequency resistance, R_{HF} :


- $R_{cell}(T, RH)$ {
- Membrane – dominate
 - Pt electrode – very small
 - GDE – very small
 - Pt electrode / GDE contact – very small
 - GDE/membrane interface – $f(RH, T)$; can be significant

➤ Non-membrane ohmic resistances (R_{cell}) must be accounted [1]

1. Cooper, *J. Electrochem. Soc.* **157**, B1731 (2010)

Comparing through-plane & in-plane conductivity of dispersion-cast Nafion[®] NR-212

$$\sigma_{\text{in-plane}} \cong \sigma_{\text{through-plane}}$$




Is the conductivity of dispersion-cast Nafion® NR-212 isotropic? ...

YES

➤ $\sigma_{IP} : \sigma_{TP} \cong 1$

- 3 temperatures
- Low to high RH


% RH	30 °C	80 °C	120 °C
20	1.30	0.87	1.20
40	1.03	0.94	1.11
60	1.10	0.95	0.92
80	1.03	0.93	0.86
90	1.00	0.90	0.86
95	1.05	0.98	0.81
\bar{X}_{1s}	1.07 _{0.08}	0.93 _{0.03}	0.99 _{0.14}



σ_{IP} = in-plane, σ_{TP} = through-plane

Through-plane & in-plane conductivity of extruded Nafion[®] 112


IP Longitudinal > IP Transverse > TP for extruded Nafion[®]




Is the conductivity of extruded Nafion® 112 isotropic? ... NO

IP Longitudinal > IP Transverse > TP

<i>Ratio of Conductivity</i>			
% RH	IP-L : TP	IP-L : IP-T	IP-T : TP
20	1.60	1.42	1.12
40	1.46	1.45	1.00
60	1.52	1.46	1.04
80	1.65	1.46	1.13
95	1.82	1.48	1.23
\bar{X}_{1s}	1.61 _{0.14}	1.46 _{0.02}	1.11 _{0.09}



- Test for equality of means
 - $H_0 : \mu_1 = \mu_2 ; H_1 : \mu_1 \neq \mu_2$
 - $\alpha = 5\%$
- Statistically significant difference in means for all RHs
 - except IP-Trans vs. TP @ 40% RH



Effective conductivity (σ_{eff}) of membrane with regions of unequal conductivity, e.g., ionomer-impregnated non-conductive porous support


$$\sigma = \frac{1}{R} \frac{L}{A}$$

In-Plane (IP)



$f \cdot Z$
 $(1-f) \cdot Z$

Through-Plane (TP)



f = fractional thickness of layer 1

In-Plane (IP)

Resistances in parallel

$$\frac{1}{R_{total}} = \frac{1}{R_1} + \frac{1}{R_2}$$

$$\sigma_{eff,IP} = f \cdot \sigma_1 + (1-f) \cdot \sigma_2$$


Through-Plane (TP)

Resistances in series


$$R_{total} = R_1 + R_2$$

$$\sigma_{eff,TP} = \frac{\sigma_1 \cdot \sigma_2}{(1-f) \cdot \sigma_1 + f \cdot \sigma_2}$$

Effective conductivity with phases of unequal conductivity



Effective conductivity with phases of unequal conductivity



$$\sigma_{eff,IP} = f \cdot \sigma_1 + (1-f) \cdot \sigma_2$$

$$\sigma_{eff,TP} = \frac{\sigma_1 \cdot \sigma_2}{(1-f) \cdot \sigma_1 + f \cdot \sigma_2}$$

f = fractional thickness of phase 1

- $\sigma_{eff, in-plane} > \sigma_{eff, through-plane}$ for supported membrane
- $\sigma_{eff, in-plane} : \sigma_{eff, through-plane}$ is a maximum for $f = 0.5$
- $\sigma_{eff, in-plane} : \sigma_{eff, through-plane} \rightarrow 1$ as $f \rightarrow 0$ or 1
- $\sigma_{eff, in-plane} : \sigma_{eff, through-plane}$ increases as $\sigma_1 : \sigma_2 \rightarrow 0$ or $\gg 1$


Comparing through-plane & in-plane conductivity of PFSA-based membrane with inert support *GORE-SELECT*[®]

- $\sigma_{eff, in-plane} > \sigma_{eff, through-plane}$ ✓
- $\sigma_{eff, in-plane} : \sigma_{eff, through-plane} \rightarrow 1$ as $f \rightarrow 0$ or 1 ✓
- ✓ Ratio is greater for thin membrane with same support thickness

Thickness $\frac{\sigma_{eff, in-plane}}{\sigma_{eff, through-plane}}$

18 μm 1.53 ± 0.16

35 μm 1.11 ± 0.10



GORE-SELECT, GORE and designs are trademarks of W. L. Gore & Associates, Inc.

Conclusions

- Methods exist that can differentiate in-plane and through-plane ionic conductivity of PEMs
- Extruded Nafion[®] 112 exhibits anisotropic conductivity
 - $\sigma_{IP-Long} > \sigma_{IP-Trans} > \sigma_{TP}$
- Dispersion-cast Nafion[®] NR-212 exhibits isotropic conductivity
 - $\sigma_{IP} \cong \sigma_{TP}$
- Membranes with inert support, *e.g.*, GORE-SELECT[®], exhibit anisotropic conductivity
 - Consistent with simple analytical treatment

Results highlight the need to consider potential for anisotropic behavior and measure appropriately

GORE-SELECT, GORE and designs are trademarks of W. L. Gore & Associates, Inc.