


# Through-Plane & In-Plane Conductivity of Polymer Electrolyte Membranes

Scribner Associates, Inc.

[www.scribner.com](http://www.scribner.com)

September 2010


# Large discrepancy in reported membrane conductivity data highlights need for accurate, robust measurement methods



# Objective: Develop an accurate & reliable test apparatus & method for *through-plane membrane resistance* & *conductivity* measurements


## ➤ Key desirable features

- ✓ Uses bare (non-catalyzed) membrane relevant thicknesses – 10 to 200  $\mu\text{m}$
- ✓ Operate over a wide range of conditions
  - 30 to  $> 120$   $^{\circ}\text{C}$
  - dry to  $> 95\%$  RH
  - 1 to 3  $\text{atm}_a$
- ✓ Rapid  $\sim 15$  min per test condition
- ✓ Robust - accurate, repeatable and reliable




$$\sigma_{\text{membrane}} = \frac{L}{R_{\text{membrane}} \cdot A} \quad [S / cm]$$

# Membrane Test System *MTS 740*



U.S. Patent No. 7,652,479


# Wet-dry gas mixing for rapid RH cycling in both directions



# Wet-dry gas mixing for controlled, rapid RH cycling


- Repeatable, reproducible and stable T, dew point & RH
  - ✓  $\pm 2\%$  from 20% to 95% RH
- Rapid RH cycling  $\rightarrow$  time-efficient testing over wide RH range
- Dew point to 120 °C, sample to 150 °C

Nominal % RH	Difference from Nominal, % RH	
	30 °C N = 15	80 °C N = 11
<b>20</b>	0.5	-0.3
<b>40</b>	0.7	-0.5
<b>60</b>	-0.2	-0.8
<b>80</b>	-1.4	-0.4
<b>90</b>	-1.5	0.3
<b>95</b>	-2.2	-0.3



# Comparing through-plane & in-plane conductivity of Nafion® NR-212

$\sigma_{\text{in-plane}} \cong \sigma_{\text{through-plane}}$  for dispersion cast Nafion®




# Is the conductivity of Nafion<sup>®</sup> isotropic? ... No consensus in published literature

$\sigma_{\parallel}$  = in-plane,  $\sigma_{\perp}$  = through-plane

$\sigma_{\parallel}:\sigma_{\perp} \cong 1$ , NR-212, this work

% RH	30 °C	80 °C	120 °C
20	1.19	0.92	1.20
40	0.95	1.02	1.15
60	1.15	1.02	0.96
80	0.99	0.96	0.87
90	0.93	0.97	0.87
95	0.93	106	0.84


- Yes, it is isotropic,  $\sigma_{\parallel}:\sigma_{\perp} \cong 1$  [Nouel; Silva]
  - ✓ This work for NR-212
- No, it is anisotropic
  - ✓  $\sigma_{\parallel}:\sigma_{\perp} = 3.6$  [Gardner]
  - ✓  $\sigma_{\parallel}:\sigma_{\perp} = 2.5 - 5$  (with pressure) [Ma]
  - ✓  $\sigma_{\parallel}:\sigma_{\perp} = 1.8 - 5$  [Casciola]
- Discrepancy due
  - ✓ Different water content ( $\lambda$ )
  - ✓ Extruded (N11X) vs. dispersion cast (NR-21X)



Gardner *et. al.*, *J. Electroanal Chem* **449** 209 (1998)  
 Ma *et. al.*, *JES* **153** A2274 (2006)  
 Casciola *et. al.*, *J. Power Sources* **162** 141 (2006)  
 Nouel, Fedkiw, *Electrochimica Acta* **43** 2381 (1998)  
 Silva, *et. al.*, *J. Power Sources* **134** 18 (2004)



# Effective conductivity ( $\sigma_{eff}$ ) of membrane with phases of unequal conductivity, e.g., ionomer-impregnated non-conductive porous support



$$\sigma_{eff, in-plane} = f \cdot \sigma_1 + (1-f) \cdot \sigma_2$$

$$\sigma_{eff, through-plane} = \frac{\sigma_1 \cdot \sigma_2}{(1-f) \cdot \sigma_1 + f \cdot \sigma_2}$$


$f$  = fractional thickness of phase 1

- $\sigma_{eff, in-plane} > \sigma_{eff, through-plane}$  for supported membrane
- $\sigma_{eff, in-plane} : \sigma_{eff, through-plane}$  is a maximum for  $f = 0.5$
- $\sigma_{eff, in-plane} : \sigma_{eff, through-plane} \rightarrow 1$  as  $f \rightarrow 0$  or  $1$
- $\sigma_{eff, in-plane} : \sigma_{eff, through-plane}$  increases as  $\sigma_1 : \sigma_2 \rightarrow 0$  or  $\gg 1$

# Comparing through-plane & in-plane conductivity ( $\sigma$ ) of PFSA-based membranes with inert support *GORE-SELECT*<sup>®</sup>

- $\sigma_{eff, in-plane} > \sigma_{eff, through-plane}$  ✓
- $\sigma_{eff, in-plane} : \sigma_{eff, through-plane} \rightarrow 1$  as  $f \rightarrow 0$  or  $1$  ✓
- ✓ Ratio is greater for thin membrane with same support thickness

Thickness	$\frac{\sigma_{eff, in-plane}}{\sigma_{eff, through-plane}}$
18 $\mu\text{m}$	1.53 $\pm$ 0.16
35 $\mu\text{m}$	1.11 $\pm$ 0.10



GORE-SELECT, GORE and designs are trademarks of W. L. Gore & Associates, Inc.



# Conclusions – Membrane Test System **MTS 740**

- ***Through-plane resistance & conductivity test system developed***
  - ✓ Bare membrane – rapid, lower cost assessment vs. MEA / fuel cell testing
  - ✓ Repeatable, accurate control of environmental conditions: cell to 150 °C, humidifier to 120 °C, dry to >95% RH
  - ✓ Robust method – repeatable and accurate
- ***Correction for non-membrane ohmic resistance contributions is important, especially for thin membranes with low resistance / high conductivity***
- ***Dispersion cast Nafion<sup>®</sup> NR-212 through-plane conductivity is the same as in-plane***
- ***Differentiate in-plane and through-plane conductivity for anisotropic material, e.g., GORE-SELECT<sup>®</sup> supported membrane***

GORE-SELECT, GORE and designs are trademarks of W. L. Gore & Associates, Inc.

Proprietary Data – Disclosure permitted only with prior written permission from Scribner Associates, Inc.


# Supporting Information

- Electrode design
- Test procedure
- Analysis Procedure
- Determination of Cell Resistance

# 4-Electrode, "Offset" Electrode Design

Side View

Dimensions in mm



Top View




U.S. Patent No. 7,652,479

# Offset electrode geometry for 4-probe conductivity measurement of thin film electrolytes

Model condition: constant 1 V between source electrodes


Source Electrode #1

100  $\mu\text{m}$



$V_{S2} = 0.5 \text{ V}$

$E, \text{ V}$



X, mm

Source Electrode #2

$V_{S1} = -0.5 \text{ V}$

Y-axis is expanded 10x relative to X-axis

U.S. Patent No. 7,652,479

# Procedure – Pre-test

- As-received membrane, stored at ambient conditions
- 32 mm x 10 mm sample
- Measure “dry” membrane thickness
  - ✓ Mean of 5 locations, 3x measurements/location
  - ✓ Low load, high accuracy gage
- Cell Assembly
  - ✓ GDE (E-LAT) cut with jig
  - ✓ Glue GDE to Pt electrode with carbon paste
  - ✓ Load membrane between GDE-prepared plattens
  - ✓ Compress ~ 2,200 kPa (325 psi) using spring loaded cell head (dial gage)






# Procedure

- Temperature series (°C): 80 → 30 → 120
- Per temperature
  - ✓ Wet-up 2 hr @ 70% RH
  - ✓ RH cycle: 70 → 20 → 90 → 95 %, 15 min step
  - ✓ Impedance sweep after 15 min
- ~ 1 day/temperature, ~6 hr
- Gas: H<sub>2</sub> or N<sub>2</sub>
- Impedance Measurement
  - ✓ 4-electrode, 4 terminal
  - ✓ Solartron 1260 FRA (standalone) / ZPlot<sup>®</sup>
  - ✓ 10 MHz – 1 Hz, 10 mV<sub>AC</sub>, 0 V<sub>DC</sub>, 10 steps/dec (~ 2 min)







Temp, °C	Total Dry Gas Flow, sccm	Pressure, kPa <sub>a</sub>
30	500	100
80	500	100
120	500	230



# Post-test Procedure – EIS Analysis



# Post-test Procedure – EIS Analysis



Element	Freedom	Value	Error	Error %
R_b	Free(+)	0.068111	0.0011616	1.7055
R1	Free(+)	12.92	0.087668	0.67854
CPE1-T	Free(+)	1.3459E-05	3.089E-07	2.2951
CPE1-P	Free(+)	0.83049	0.0017289	0.20818

# Through-plane Resistance & Conductivity

- Through-plane resistance includes ***non-membrane ohmic contributions,  $R_{cell}$***

$$\text{Cell Resistance, } R_{cell} = R_{\Omega, electrode} + R_{\Omega, contact} + R_{\Omega, interface}$$

- Typically work in area specific resistance, ***ASR***

$$ASR_{uncorrected} = R_{HF} \cdot A_{effective} \quad [\Omega - cm^2]$$

- Accounting for the cell ***ASR(T, RH)*** gives the membrane resistance

$$ASR_{membrane}(T, RH) = ASR_{uncorrected}(T, RH) - ASR_{cell}(T, RH) \quad [\Omega - cm^2]$$


- ✓ Note that all are a ***f(T, RH)***

**The challenge: need  $ASR_{cell}(T, RH)$**

# Determine $ASR_{cell}$ by extrapolating linear regression of ASR vs. thickness to $L = 0$ . Do this for each T, RH

- $R^2 = 0.95 - 0.99$
- Similar  $ASR_{cell}$  for 2 sets of PFSA membranes
  - ✓ Nafion® N1XX (4 thicknesses)
  - ✓ Supported PFSA membrane (3 thicknesses)

$ASR_{cell}$




**Key Assumptions:**

1. Intrinsic through-plane conductivity is not a function of  $L$
2. Cell resistance is constant from build-to-build

# Cell resistance increases at low RH ... dominated by interfacial resistance

➤ Increasing  $R_{\text{cell}}$  with decreasing RH also reported by W.L. Gore & Associates <sup>1</sup>



1. Johnson & Liu, "Ionic Conductivity of Perfluorosulfonic Acid Membranes as a Function of Temperature, Humidity and Equivalent Weight" ECS PV 2002-5, 132 (2002)

# Ratio of cell to membrane resistance highlights importance of correcting for non-membrane ohmic contributions

- At low RH, ratio is small relative to the membrane resistance
  - ✓  $R_{\text{cell}} : R_{\text{membrane}} \sim 0.1 - 0.2$
- At high RH, the cell resistance can be significant relative to the membrane resistance, especially for thin membranes

$$\frac{\text{Cell } R}{\text{Membrane } R}$$

